



Outdoor Technical Emergency Care

Better Solutions. Better Care.



Aer Cervical Collar

Packable, lightweight and durable cervical collar.

User friendly design for easy and efficient placement.

KEY FEATURES

Chin strap provides intuitive placing and comfort of collar.

Easy and secure **latching** allows for varying sizes.

Air battens create structure while maintaining critical blood flow to the brain.

Durable 75d Ripstop Nylon is strong for all situations and environments.

Packable size allows for easy carry and storage (2.4 oz.).



CONTACT US TODAY

801.810.0048 | admin@otecdesign.com | otecdesign.com

